

# ABI-400CL

Automated ABI Testing

CPT 93922

# Newman Medical

simpleABI Systems from Newman Medical use the power of computers to produce the reimbursement reports you need quickly, easily, and accurately. The easy to use PC-based systems' minimal learning curve makes it easy to begin performing exams. Automatic index calculations make the simpleABI system make exams as easy as taking blood pressures. **Print to any standard office printer or save and load directly to EMR.** The diagnosis is made using the gold standard Doppler and simple pulse volume recording (PVR) waveforms. The ABI-400CL system meets the new dual ankle requirements for CPT code 93922.

## ABI-400CL Features:

- **Cuff-Link™ automated cuff select, inflate, calculate**
- PC-based ABI system
  - Easy staff training - simple to learn and use
  - Point & Click reporting right on the computer
- Uses "Gold Standard" 8MHz Doppler
- Easy "Click & Capture" PVR waveforms
- Customizable reports with practice information and logo
- **Print to any standard office printer**
- **Attach multiple formats directly to EMR records**
- Full 2-year product warranty

## ABI-400CL Kit Contains:

- Cuff-Link™-4
- Single Level ABI & TBI Testing
- Sensitive 8MHz Doppler
- PPG Probe
- 5 cuffs (2x10cm, 2x12cm, 1x1.9cm)



Cuff-Link – automatic cuff selection



Cuff-Link Remote - 6 button simplicity!

## ABI-400CL Options:

- Computer
- Roll Stand
- Printer
- Mobile Carry Bag



Newman Medical  
5350 Vivian Street | Unit C | Arvada | CO | 80002  
Toll-free 800-267-5549 | Fax-888-202-3627  
[www.newman-medical.com](http://www.newman-medical.com)

NM0004C