

# ABI-500CL

Automated Single & Multi-level Testing  
CPT 93922 & 93923

# Newman Medical

simpleABI Systems from Newman Medical use the power of computers to produce the reimbursement reports you need quickly, easily, and accurately. The PC-based systems' minimal learning curve makes it easy to begin performing exams. Automatic index calculations make the simpleABI system make exams as easy as taking blood pressures. **Print to any standard office printer or save and load directly to EMR.** The diagnosis is made using the gold standard Doppler and simple pulse volume recording (PVR) waveforms. The ABI-500CL system meets all requirements, including dual ankle, for CPT code 93922 & 93923.

## ABI -500CL Features:

- **Cuff-Link™** automated cuff select, inflate, calculate
- PC-based ABI system
  - Easy staff training - simple to learn and use
- Uses "Gold Standard" 8MHz Doppler
- Easy "Click & Capture" PVR waveforms
- Customize reports with practice information and logo
- **Print to any standard office printer**
- **Attach multiple formats directly to EMR records**
- Full 2-year product warranty

## ABI -500CL Kit Contains:

- Cuff-Link™ -8
- Single Level ABI & TBI Testing
- **Multi-level Segmental Testing**
- Sensitive 8MHz Doppler
- PPG Probe
- 9 cuffs (4x10cm, 4x12cm, 1x1.9cm)



## ABI -500CL Options:

- Computer
- Roll Stand
- Printer
- Mobile Carry Bag



Newman Medical  
5350 Vivian Street | Unit C | Arvada | CO | 80002  
Toll-free 800-267-5549 | Fax-888-202-3627  
[www.newman-medical.com](http://www.newman-medical.com)

NM0004C