

# 300 SERIES

## *Procedure Chair*

- Powered height, backrest, tilting & leg rest
- Programmable memory positions
- 550 lb patient capacity



Shown in Pewter fabric with Grey Enclosure, Stirrups, T-Rail, and Adjustable Head Rest & Arm Rest accessories.



ISO 13485:2003  
FM 584901

made in the USA  
with US & imported parts



**OAKWORKS<sup>®</sup>**  
**MEDICAL**  
*Position for Success*

# 300 SERIES *Procedure Chair*

The 300 Series Procedure Chair gives you the accessibility, positioning and patient comfort for a wide range of procedures in a compact and aesthetically pleasing package. This power chair features all electric movements with memory presets and a wide array of options and accessories.



## 300 SERIES PROCEDURE CHAIR #73076-T-E\* 3" Comfort #73408-T-E\* 4" Comfort

FEATURES	
3 Section Top	27" (69cm) wide x 73" (185cm) long with Head Rest [80" (203cm) long with Adjustable Head Rest extended]
Back rest	27" (69cm) wide x 26" (66cm) long
Seat	27" (69cm) wide x 23" (58cm) long
Leg rest	23" (58cm) wide x 13" (33cm) long
Motion 1	Electric Height Range: 19" - 35" (48 - 90cm)
Motion 2	Electric Trendelenburg/Reverse Tilt: +15°/-5°
Motion 3	Electric Backrest: 0-80°
Motion 4	Electric Leg Rest: 0-90°
Hand Control	All motions and memory positions
Weight Capacity	550 lb (250kg)
Product Weight	300 lb (136kg)
Padding	3" (8cm) or 4" (10cm) Comfort Foam™
Fabric	Premium, PVC free TerraTouch™
Enclosure	Choose from Black, Grey, Beige, or CoolWhite
Integrated Power Outlet	USA model
Paper Roll Holder	Accommodates up to 18" (46cm) rolls (rolls not included)
Warranty	3 years parts, 2 years labor
Safety Listings	FDA listed, CE marked, cETLus approved, CB Certificate

OPTIONS <i>must be ordered with table.</i>		#
Foot Control	multi-function	73074
Non-USA electric system	please contact us for order details	73134
Procedure Tray	12.25" (6cm) x 10.25" (26cm) and holds 2 quarts (1.9 liters) Available only with 4" Comfort Foam™ Reduces foam thickness in seat section to 2"	73042
Stirrups (Retractable Foot Rests)	Pull out from under seat section, fully adjustable	73175
T-rail, Fixed	Available for factory mount on one or both sides of the seat. 14" (36cm) long	73043
ACCESSORIES <i>order any time.</i>		#
Arm Rests	Removable arm rests rotate horizontally and always stay level. 7.5" (19cm) wide x 23" (58cm) long	73064-T*
Adjustable Head Rest Platform	Adjustable angle & height, removable	71899
Procedure Chair Face Rest Crescent	For use when prone positioning is required. Use with Adjustable Head Rest Platform	69876-T*
Head Rest Pad	10" (25cm) long x 10" (25cm) wide Use with Adjustable Head Rest Platform	71976-T*
Large Head Rest	adjustable in length, removable 11.5" (29cm) long x 20" (51cm) wide	72717-T*
T-rail Adapter	Attaches to chair at multiple locations for use with USA format T-rail attachable accessories.	61902

\*Fabric code (T) and Enclosure color code (E) must be specified when ordering.

## ADJUSTABLE HEAD REST PLATFORM & PAD

#71899 #71976-T\*

This removable head rest is adjustable in angle, height, and distance from the table.

10" (25cm) long x 10" (25cm) wide



ACCESSORY

## T-RAIL, FIXED #73043

Factory mounted 14" (36cm) long t-rail located on either side of the seat.



OPTIONAL

## LARGE HEAD REST #72717-T\*

This accessory is available by itself without the adjustable head rest platform and the head rest pad. It is removable and is adjustable in length. 11.5" (29cm) long x 20" (51cm) wide



ACCESSORY

## 300 SERIES ARM REST #73064-T\*

These incredibly versatile, removable arm rests rotate horizontally and always stay level. 7.5" (19cm) wide at widest point x 23" (58cm) long



ACCESSORY

## FOOT CONTROL #73074

Add a foot control to adjust height, backrest, leg rest, and trendelenburg.



OPTIONAL

## T-RAIL ADAPTER #61902

Multiple locations allow you to choose the best spot for this easily installed & removed accessory capable of handling up to 150 lb (68kg). Can be used with any American standard accessory. 8" (20cm) long



ACCESSORY

## STIRRUPS #73175

Adjust to a wide range of angles & lengths to accommodate any patient or procedure. Retractable design stores inside the table.



OPTIONAL

## PROCEDURE TRAY #73042

Retractable, removable stainless steel.

12.25" (6cm) x 10.25" (26cm) and holds 2 quarts (1.9 liters)



OPTIONAL

## ENCLOSURE COLORS

Choice of 4 acrylic colors to match your decor. \*Color code must be specified.



## PADDING



**3" Comfort Foam™**  
 More comfort and support



**4" Comfort Foam™**  
 Superior comfort and support

## TERRATOUGH™ FABRIC

Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables and passes testing for Cytotoxicity, Skin Irritation and Sensitization as well as CA TB117 Fire Retardance testing. For disinfecting purposes you may use Protex, MadaCide, MadaWipe, Accell TB, Virox® (Canada), PDI Duper Sani-Cloth Germicidal wipe or a diluted 1/10 bleach solution.



We strive to ensure that the colors displayed on our product imagery (online and in print) are as accurate as possible. However, due to the limitations of printers, scanners, computer monitors, and other equipment, the colors may not be accurate.

CONTACT:



*View our full  
 line of Medical  
 products #72703*



Scan this code  
 with your  
 smart phone